

ZytoLight® SPEC Williams-Beuren Dual Color Probe



Background

The ZytoLight® SPEC Williams-Beuren Dual Color Probe is designed to detect deletions affecting the chromosomal region 7q11.23 harboring the ELN (elastin, a.k.a. WBS) gene.

The Williams-Beuren syndrome (WBS) is a genetic disorder caused by a hemizygous contiguous gene deletion on chromosome 7q11.23. The estimated prevalence of the disease ranges between 1/7,500 and 1/20,000 newborns.

The WBS deletion region (~1.5-1.8 Mb) consists of a single copy gene region containing app. 28 genes, including the ELN gene that is flanked by repetitive sequences known as low-copy repeats (LCRs).

The deletions arise as a consequence of misalignment of these repetitive sequences during meiosis and a following unequal crossing over due to high similarity of LCRs. Usually, WBS occurs sporadically, but some parents of WBS patients were shown to carry a paracentric inversion of the WBS locus. Presence of this inversion predisposes to chromosomal mispairing in meiosis.

WBS patients clinically display a characteristic pattern of symptoms including vascular stenosis, weakness of connective tissue, a typical face, short stature, over-friendliness, and mental retardation. FISH analysis can be performed to confirm WBS diagnosis in patients with vascular stenosis together with mental retardation.


References

- Bayés M, et al. (2003) Am J Hum Genet 73: 131-51.
- Beuren AJ, et al. (1964) Am J Cardiol 13: 471-83.
- Schubert C (2009) Cell Mol Life Sci 66: 1178-97.
- Sugayama SM, et al. (2003) Arq Bras Cardiol 81: 462-73.
- Williams JC, et al. (1961) Circulation 24: 1311-8.


Probe Description

The ZytoLight® SPEC Williams-Beuren Dual Color Probe is composed of:

- ZyOrange (excitation 547 nm/emission 572 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 7q11.23** (chr7:73,408,390-73,899,599) harboring the ELN gene region.
- ZyGreen (excitation 503 nm/emission 528 nm) labeled polynucleotides (~4.5 ng/µl), which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
- Formamide based hybridization buffer




Ideogram of chromosome 7 indicating the hybridization locations.




SPEC ELN Probe map (not to scale).

Results

In a normal interphase nucleus, two orange and two green signals are expected. In a cell with deletion of the ELN gene locus, a reduced number of orange signals will be observed.



SPEC Williams-Beuren Dual Color Probe hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus and to metaphase chromosomes of a normal cell.



Lymphocytes and metaphase chromosomes from a Williams-Beuren syndrome case showing an ELN deletion as indicated by the loss of one orange signal.

Prod. No.	Product	Label	Tests* (Volume)
Z-2302-50	ZytoLight SPEC Williams-Beuren Dual Color Probe	●/●	5 (50 µl)
Related Products			
Z-2099-20	ZytoLight FISH-Cytology Implementation Kit Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl ₂ , 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTest-Solution, 0.8 ml		20

* Using 10 µl probe solution per test. labeled products are only available in certain countries. All other countries research use only! Please contact your local dealer for more information.

**According to Human Genome Assembly GRCh37/hg19