



**ZytoFast**  
**EBV Probe**  
(Digoxigenin-labeled)

REF T-1114-400  40 (0.4 ml)

For the qualitative detection of  
human Epstein-Barr virus (EBV) EBER RNA  
by chromogenic *in situ* hybridization (CISH)



In vitro diagnostic medical device  
according to EU directive 98/79/EC

### 1. Intended use

The ZytoFast EBV Probe (PF29) is intended to be used for the qualitative detection of human Epstein-Barr virus (EBV) EBER RNA in formalin-fixed, paraffin-embedded specimens, such as lymphomas, by chromogenic *in situ* hybridization (CISH). The probe is intended to be used in combination with the ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40).

Interpretation of the results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

### 2. Clinical relevance

EBV (a.k.a. human herpesvirus-4, HHV-4) is a member of the gamma-herpesvirus group and one of the most common viruses in humans. Transmission of EBV requires close, intimate contact with a person excreting the virus in its saliva. EBV has two major target tissues *in vivo*, B lymphocytes and squamous pharyngeal epithelium. Infection of B lymphocytes with EBV results in persistent latent infection, immortalization of the cells, and perpetual proliferation. EBV, the first virus to be identified as an oncovirus, is the etiological agent of infectious mononucleosis and has been implicated in the pathogenesis of an increasing number of human malignancies such as Burkitt's lymphoma, nasopharyngeal carcinoma, and polyclonal lymphomas in immuno-compromised individuals. CISH-based diagnosis of EBV infection has the advantage over other methods in that it permits unequivocal localization of EBV genomes in cells and thereby obviates the risk of false positive results due to laboratory or clinical contamination.

### 3. Test principle

The chromogenic *in situ* hybridization (CISH) technique allows the detection and visualization of specific nucleic acid sequences in cell preparations. Hapten-labeled nucleotide fragments, so called CISH probes, and their complementary target sequences in the preparations are co-denatured and subsequently allowed to anneal during hybridization. Afterwards, unbound and unbound probe fragments are removed by stringency washing steps. Duplex formation of the labeled probe can be visualized using primary (unmarked) antibodies, which are detected by secondary polymerized enzyme-conjugated antibodies. The enzymatic reaction with chromogenic substrates subsequently leads to the formation of colored precipitates. After counterstaining the nucleus with a nuclear dye, hybridized probe fragments are visualized by light microscopy.

### 4. Reagents provided

The ZytoFast EBV Probe (PF29) is composed of:

- Digoxigenin-labeled oligonucleotides (~ 0.2 ng/μl), which target mRNA sequences encoding EBER-1 and EBER-2 regions.

The ZytoFast EBV Probe is available in one size:

- T-1114-400: 0.4 ml (40 reactions of 10 μl each)

### 5. Materials required but not provided

- ZytoFast PLUS CISH Implementation Kit HRP-DAB (Prod. No. T-1063-40) or
- Positive and negative control specimens
- Microscope slides, positively charged
- Water bath (55°C, 98°C)
- Hybridizer or hot plate
- Hybridizer or humidity chamber in hybridization oven
- Adjustable calibrated pipettes (10 μl, 100 μl, 1000 μl)
- Staining jars or baths
- Timer
- Calibrated thermometer
- Ethanol or reagent alcohol
- Xylene
- Methanol 100%
- Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) 30%
- Deionized or distilled water
- Coverslips (22 mm x 22 mm, 24 mm x 32 mm)
- Rubber cement, e.g., Fixogum Rubber Cement (Prod. No. E-4005-50/-125) or similar
- Adequately maintained light microscope (100-200x)

### 6. Storage and handling

Store at 2-8°C in an upright position. Return to storage conditions immediately after use. Do not use reagents beyond expiry date indicated on the label. The product is stable until expiry date indicated on the label when handled accordingly.

### 7. Warnings and precautions

- Read the instructions for use prior to use!
- Do not use the reagents after the expiry date has been reached!
- This product contains substances (in low concentrations and volumes) that are harmful to health and potentially infectious. Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- If reagents come into contact with skin, rinse skin immediately with copious amounts of water!
- A material safety data sheet is available on our homepage ([www.zytovision.com](http://www.zytovision.com)).
- Do not reuse reagents!

- Avoid any cross-contamination and micro-bacterial contamination of the reagents!
- The specimens must not be allowed to dry during the hybridization and washing steps!

#### Hazard and precautionary statements for PF29:

This probe is not classified as hazardous according to Regulation (EC) No. 1272/2008.

### 8. Limitations

- For *in vitro* diagnostic use.
- For professional use only.
- The clinical interpretation of any positive staining, or its absence, must be done within the context of clinical history, morphology, other histopathological criteria as well as other diagnostic tests. It is the responsibility of a qualified pathologist to be familiar with the CISH probes, reagents, diagnostic panels, and methods used to produce the stained preparation. Staining must be performed in a certified, licensed laboratory under the supervision of a pathologist who is responsible for reviewing the stained slides and assuring the adequacy of positive and negative controls.
- Specimen staining, especially signal intensity and background staining, is dependent on the handling and processing of the specimen prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other specimens or fluids may produce artefacts or false results. Inconsistent results may result from variations in fixation and embedding methods, as well as from inherent irregularities within the specimen.
- The probe should be used only for detecting target sequences described in 4. "Reagents provided".
- The performance was validated using the procedures described in these instructions for use. Modifications to these procedures might alter the performance and have to be validated by the user.

### 9. Interfering substances

The following fixatives are incompatible with ISH:

- Bouin's fixative
- B5 fixative
- Acidic fixatives (e.g., picric acid)
- Zenker's fixative
- Alcohols (when used alone)
- Mercuric chloride
- Formaldehyde/zinc fixative
- Hollande's fixative
- Non-buffered formalin

### 10. Preparation of specimens

Recommendations:

- Avoid cross-contamination of samples in any step of preparation as this may lead to erroneous results.
- Fixation in 10% neutrally buffered formalin for 24 h at room temperature (18-25°C).
- Sample size  $\leq 0.5 \text{ cm}^3$ .
- Use premium quality paraffin.
- Embedding should be carried out at temperatures lower than 65°C.
- Prepare 3-5  $\mu\text{m}$  microtome sections.
- Use positively charged microscope slides.
- Fix tissue sections for 2-16 h at 50-60°C.

### 11. Preparatory treatment of the device

The product is ready-to-use. No reconstitution, mixing, or dilution is required. Bring probe to hybridization temperature (55°C) and mix thoroughly before use.

## 12. Assay procedure

#### Specimen pretreatment

Perform specimen pretreatment (e.g., dewaxing, proteolysis) according to the instructions for use of the ZytoFast PLUS CISH Implementation Kits.

#### Denaturation and hybridization

1. Pipette 10  $\mu\text{l}$  of the probe onto each pretreated specimen.
2. Cover specimens with a 22 mm x 22 mm coverslip (avoid trapped bubbles) and seal the coverslip.

*We recommend using rubber cement (e.g., Fixogum) for sealing.*

3. Place slides on a hot plate or hybridizer and denature specimens for 5 min at 75°C.
4. Transfer slides to a humidity chamber and hybridize for 1 h at 55°C (e.g., in a hybridization oven).

*It is essential that specimens do not dry out during the hybridization step.*

#### Post-hybridization and detection

Perform post-hybridization processing (washing, detection, counter-staining, mounting, microscopy) according to the instructions for use of the respective ZytoFast PLUS CISH Implementation Kit.

## 13. Interpretation of results

Using the ZytoFast PLUS CISH Implementation Kit, hybridized Digoxigenin-labeled oligonucleotides appear as brown pattern when detected by horseradish peroxidase (HRP) and DAB.

A positive reactivity for Epstein-Barr virus (EBV) EBER RNA in the target cells is indicated by a distinctly stained nucleus.

#### Please note:

- Visualization of signals should be performed at least at 100-fold magnification resulting in easily visible signals.
- Do not evaluate areas of necrosis, overlapping nuclei, over-digested nuclei, and nuclei with weak signal intensity.
- A negative or unspecific result can be caused by multiple factors (see chapter 17 "Troubleshooting").
- In order to correctly interpret the results, the user must validate this product prior to use in diagnostic procedures according to national and/or international guidelines.

## 14. Recommended quality control procedures

In order to monitor correct performance of processed specimens and test reagents, each assay should be accompanied by internal and external controls. If internal and/or external controls fail to demonstrate appropriate staining, results with patient specimens must be considered invalid.

**Internal controls:** In unclear cases, RNA Control Probes should be used for further clarification.

**External controls:** Validated positive and negative control specimens.

## 15. Performance characteristics

In order to assess the analytical performance of the ZytoFast human EBV Probe, chromogenic *in situ* hybridizations were performed on formalin-fixed, paraffin-embedded (FFPE) specimens of three different cell lines. (Raji, NC-NC, and JVM-2) and a cell control array (MB-CC VIR).

**Analytical sensitivity:** The overall analytical sensitivity of the ZytoFast EBV Probe assessed by CISH on EBV-positive cell lines was calculated to be 100.0 %.

**Analytical specificity:** The overall analytical specificity of the ZytoFast EBV Probe assessed by CISH on EBV-negative specimens was calculated to be 100.0 %.

## 16. Disposal

The disposal of reagents must be carried out in accordance with local regulations.

## 17. Troubleshooting

Any deviation from the operating instructions can lead to inferior staining results or to no staining at all.

**Weak signals or no signals at all**

| Possible cause   | Action  |
|--|---|
| Cell or tissue sample has not been properly fixed  | Optimize fixing time and fixative   |
| Proteolysis, hybridization, denaturation, stringency wash or antibody-incubation temperature not correct | Check temperature of all technical devices used, using a calibrated thermometer. Use always the same number of slides in solutions with critical temperature  |
| Proteolytic pretreatment not carried out properly  | Depending on multiple factors, e.g., nature and duration of fixing, thickness of sections, and nature of tissue/cells, different incubation times may be required. Ascertain the optimum time for pepsin incubation in pre-tests  |
| Hybridization time too short   | Hybridize for at least 1 h; extend hybridization time if necessary  |
| Too low concentrated Wash Buffer   | Check concentration of Wash Buffer  |
| Old dehydration solutions  | Prepare fresh dehydration solutions   |
| Probe evaporation  | When using a hybridizer, the use of the wet stripes/water filled tanks is mandatory. When using a hybridization oven, the use of a humidity chamber is required. In addition, the coverslip should be sealed completely, e.g., with Fixogum, to prevent drying-out of the sample during hybridization |
| Insufficient preparation of chromogenic substrate  | Instead of preparing the color substrates by dropping, use a pipette  |
| Incubation temperature for color substrates not correct  | Check temperature of all technical devices used, using a calibrated thermometer   |
| Counterstaining time too long  | The counterstaining time depends on the nature of the specimen and should be optimized accordingly. Avoid dark counterstaining, because it may obscure positive staining signals  |
| No target sequences available  | Use positive control probes to verify pepsin incubation time. Use verified positive tissue to confirm test performance  |
| Bluing of counterstain not carried out properly  | Use cold running tap water for bluing; do not use warm or hot water or bluing reagents  |

**Signals too strong**

| Possible cause                                | Action   |
|---|--|
| Proteolytic pretreatment carried out too long | Depending on multiple factors, e.g., nature and duration of fixing, thickness of sections, and nature of tissue/cells, different incubation times may be required. Ascertain the optimum time for pepsin incubation in pre-tests |
| Substrate reaction is too intense             | Shorten substrate incubation time; do not heat substrate solution above the temperature given in the instructions for use  |

**Signals fade or merge**

| Possible cause                                | Action   |
|---|--|
| An unsuitable mounting solution has been used | Use only the mounting solution provided with the kit or recommended by the instructions for use. Use solutions free of any impurities; do not use coverslip tape |

**Uneven or in some parts only very light staining**

| Possible cause   | Action  |
|--|---|
| Incomplete dewaxing  | Use fresh solutions; check duration of dewaxing                         |
| Reagent volume too small                                   | Ensure that the reagent volume is large enough to cover the tissue area |
| Air bubbles caught before hybridization or during mounting | Avoid air bubbles   |

**Inconsistent results**

| Possible cause  | Action   |
|---|--|
| Insufficient drying before probe application  | Extend air-drying  |
| Too much water/wash buffer on tissue prior to application of pepsin, antibodies and/or color substrates | Ensure that excess liquid is removed from tissue section by blotting or shaking it off the slide. Small amounts of residual water/wash buffer do not interfere with the test |
| Variations in tissue fixation and embedding methods   | Optimize fixation and embedding methods  |
| Variations in tissue section thickness  | Optimize sectioning  |

**Tissue morphology degraded**

| Possible cause                                    | Action                            |
|---|-----------------------------------|
| Cell or tissue sample has not been properly fixed | Optimize fixing time and fixative |
| Proteolytic pretreatment not carried out properly | Optimize pepsin incubation time   |

### Noisy background

| Possible cause  | Action   |
|---|--|
| Stringency wash temperature not correct                   | Check temperature of the technical devices used, using a calibrated thermometer. Use always the same number of slides in the jar. We recommend not to use more than eight slides per jar for heat incubation steps |
| Slides not thoroughly rinsed                              | Use fresh and sufficient wash buffer and deionized or distilled water where indicated  |
| Sections dried out any time during or after hybridization | Avoid sections being dried out; use humidity chamber; seal coverslip properly  |
| Prolonged substrate incubation time                       | Shorten substrate incubation time  |
| Incomplete dewaxing                                       | Use fresh solutions; check duration of dewaxing  |
| Proteolytic pretreatment too strong                       | Optimize pepsin incubation time  |
| Slides cooled to room temperature before hybridization    | Transfer the slides quickly to hybridization temperature   |
| Tissue-antibody interaction                               | Use negative control probes to ascertain tissue-specific background staining   |

### Overlapping signals

| Possible cause                             | Action                                       |
|--|--|
| Inappropriate thickness of tissue sections | Prepare 3-5 $\mu\text{m}$ microtome sections |

### Specimen floats off the slide

| Possible cause                      | Action                                      |
|-------------------------------------|---|
| Unsuitable slide coating            | Use appropriate (positively charged) slides |
| Proteolytic pretreatment too strong | Shorten pepsin incubation time              |

## 18. Literature

- Greifenegger N, et al. (1998) *J Virol* 72: 9323-8.
- Khan G, et al. (1992) *J Clin Pathol* 45: 616-20.
- Khanna R, et al. (1995) *Microbiol Rev* 59: 387-405.
- Murphy JK, et al. (1990) *J Clin Pathol* 43: 220-3.
- Rosa MD, et al. (1981) *Mol Cell Biol* 1: 785-96.
- Thorley-Lawson DA (2001) *Nat Rev Immunol* 1: 75-82.
- Wilkinson DG: *In Situ Hybridization, A Practical Approach*, Oxford University Press (1992), ISBN 0 19 963327 4.

Our experts are available to answer your questions.  
Please contact [help@zytovision.com](mailto:help@zytovision.com)



ZytoVision GmbH  
Fischkai 1  
27572 Bremerhaven/ Germany  
Phone: +49 471 4832-300  
Fax: +49 471 4832-509  
[www.zytovision.com](http://www.zytovision.com)  
Email: [info@zytovision.com](mailto:info@zytovision.com)

### Trademarks:

ZytoVision® and ZytoFast® are trademarks of ZytoVision GmbH.