



## ZytoFast PLUS

### CISH Implementation Kit AP-Permanent Red

REF T-1151-40

40

For the qualitative detection of Digoxigenin-labeled  
ZytoFast Probes by chromogenic *in situ* hybridization  
(CISH)

For research use only.  
Not for use in diagnostic procedures.

#### 1. Intended use

The ZytoFast PLUS CISH Implementation Kit AP-Permanent Red is intended to be used in combination with Digoxigenin-labeled ZytoFast Probes on formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

#### 2. Test principle

The chromogenic *in situ* hybridization (CISH) technique allows the detection and visualization of specific nucleic acid sequences in cell preparations. Hapten-labeled nucleotide fragments, so called CISH probes, and their complementary target sequences in the preparations are co-denatured and subsequently allowed to anneal during hybridization. Afterwards, unspecific and unbound probe fragments are removed by stringency washing steps. Duplex formation of the labeled probe can be visualized using primary (unmarked) antibodies, which are detected by secondary polymerized enzyme-conjugated antibodies. The enzymatic reaction with chromogenic substrates leads to the formation of colored precipitates. After counterstaining the nucleus with a nuclear dye, hybridized probe fragments are visualized by light microscopy.

#### 3. Reagents provided

The ZytoFast PLUS CISH Implementation Kit AP-Permanent Red is available in one size and is composed of:

| Code | Component                              | Quantity<br>40 | Container                  |
|------|--|----------------|----------------------------|
| PT2  | <u>Heat Pretreatment Solution EDTA</u> | 500 ml         | Screw-cap bottle (large)   |
| ES1  | <u>Pepsin Solution</u>                 | 4 ml           | Dropper bottle, white cap  |
| WB5  | <u>20x Wash Buffer TBS</u>             | 4x 50 ml       | Screw-cap bottle           |
| AB11 | <u>Rabbit-Anti-Dig</u>                 | 4 ml           | Dropper bottle, grey cap   |
| AB12 | <u>Anti-Rabbit-AP-Polymer</u>          | 4 ml           | Dropper bottle, green cap  |
| SB8a | <u>Permanent Red Solution A</u>        | 0.25 ml        | Reaction vessel, brown lid |
| SB8b | <u>Permanent Red Solution B</u>        | 15 ml          | Dropper bottle, red cap    |
| CS1  | <u>Mayer's Hematoxylin Solution</u>    | 20 ml          | Screw-cap bottle, black    |
| MT4  | <u>Mounting Solution (alcoholic)</u>   | 4 ml           | Glass bottle, brown        |
|      | Instructions for use                   | 1              |                            |

**T-1151-40 (40 tests):** Components **ES1**, **AB11**, **AB12**, **SB8a-b**, **CS1**, and **MT4** are sufficient for 40 reactions. Component **PT2** is sufficient for 7 staining jars of 70 ml each. Component **WB5** is sufficient for 57 staining jars of 70 ml each.

#### 4. Materials required but not provided

- Digoxigenin-labeled ZytoFast CISH Probe
- Positive and negative control tissue
- Microscope slides, positively charged
- Water bath (55°C, 98°C)
- Hybridizer or hot plate
- Hybridizer or humidity chamber in hybridization oven
- Adjustable calibrated pipettes (10 µl, 1000 µl)
- Staining jars or baths
- Timer
- Calibrated thermometer
- Ethanol or reagent alcohol
- Xylene
- Deionized or distilled water
- Coverslips (22 mm x 22 mm, 24 mm x 32 mm)
- Rubber cement, e.g., Fixogum Rubber Cement (Prod. No. E-4005-50/-125) or similar
- Adequately maintained light microscope (100-200x)

#### 5. Storage and handling

Store at 2-8°C in an upright position. Return to storage conditions immediately after use. Do not use reagents beyond expiry date indicated on the label. The product is stable until expiry date indicated on the label when handled accordingly.

#### 6. Warnings and precautions

- Read the instructions for use prior to use!
- Do not use the reagents after the expiry date has been reached!
- This product contains substances (in low concentrations and volumes) that are harmful to health and potentially infectious. Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- If reagents come into contact with skin, rinse skin immediately with copious amounts of water!
- A material safety data sheet is available on our homepage ([www.zytovision.com](http://www.zytovision.com)).
- Do not reuse reagents.
- Avoid any cross-contamination and micro-bacterial contamination of the reagents!
- The specimens must not be allowed to dry during the hybridization and washing steps!

#### Special labeling of ES1:

EUH208 Contains Pepsin A. May produce an allergic reaction.  
EUH210 Safety data sheet available on request.

#### Special labeling of CS1:

EUH210 Safety data sheet available on request.

**Hazard and precautionary statements for AB11, AB12, PT2 and WB5**

The hazard-determining component is a mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1).

**Warning**

|           |  |
|-----------|--|
| H317      | May cause an allergic skin reaction.                                       |
| P261      | Avoid breathing dust/fume/gas/mist/vapours/spray.                          |
| P272      | Contaminated work clothing should not be allowed out of the workplace.     |
| P280      | Wear protective gloves/protective clothing/eye protection/face protection. |
| P302+P352 | IF ON SKIN: Wash with plenty of water.                                     |
| P333+P313 | IF skin irritation or rash occurs: Get medical advice/attention.           |
| P362+P364 | Take off contaminated clothing and wash it before reuse.                   |

**Hazard and precautionary statements for MT4**

The hazard-determining component is xylene.

**Warning**

|                |  |
|----------------|--|
| H226           | Flammable liquid and vapour.   |
| H312+H332      | Harmful in contact with skin or if inhaled.  |
| H315           | Causes skin irritation.  |
| H319           | Causes serious eye irritation.   |
| H335           | May cause respiratory irritation.  |
| H373           | May cause damage to organs through prolonged or repeated exposure.   |
| P210           | Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.                                   |
| P260           | Do not breathe dust/fume/gas/mist/vapours/spray.   |
| P280           | Wear protective gloves/protective clothing/eye protection/face protection.   |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313      | IF eye irritation persists: Get medical advice/attention.  |
| P403+P235      | Store in a well-ventilated place. Keep cool.   |
| EUH208         | Contains methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate. May produce an allergic reaction.         |

**Hazard and precautionary statements for SB8a**

The hazard-determining component is hydrogen chloride.

**Danger**

|                |  |
|----------------|--|
| H314           | Causes severe skin burns and eye damage.   |
| H332           | Harmful if inhaled.  |
| P260           | Do not breathe dust/fume/gas/mist/vapours/spray.   |
| P280           | Wear protective gloves/ protective clothing/ eye protection/ face protection.  |
| P301+P330+P331 | IF SWALLOWED: rinse mouth. Do NOT induce vomiting.   |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.                           |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P310           | Immediately call a POISON CENTER/doctor  |

**7. Limitations**

- For research use only.
- For professional use only.
- Specimen staining, especially signal intensity and background staining, is dependent on the handling and processing of the specimen prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other specimens or fluids may produce artefacts or false results. Inconsistent results may result from variations in fixation and embedding methods, as well as from inherent irregularities within the specimen.
- Probes should be used only for detecting target sequences described in the instructions for use of the respective probe.
- The performance was validated using the procedures described in these instructions for use. Modifications to these procedures might alter the performance and have to be validated by the user.

**8. Interfering substances**

The following fixatives are incompatible with ISH:

- Bouin's fixative
- B5 fixative
- Acidic fixatives (e.g., picric acid)
- Zenker's fixative
- Alcohols (when used alone)
- Mercuric chloride
- Formaldehyde/zinc fixative
- Hollande's fixative
- Non-buffered formalin

**9. Preparation of specimens**

Recommendations:

- Avoid cross-contamination of samples as this may lead to erroneous results.
- Fixation in 10% neutrally buffered formalin for 24 h at room temperature (RT, 18°C-25°C).
- Sample size  $\leq 0.5 \text{ cm}^3$ .
- Use premium quality paraffin.
- Embedding should be carried out at temperatures lower than 65°C.
- Prepare 3-5  $\mu\text{m}$  microtome sections.
- Use positively charged microscope slides.
- Fix tissue sections for 2-16 h at 50-60°C.

**10. Preparatory treatment of the device**

20x Wash Buffer TBS (WB5) is to be prepared according to the instructions in 12. "Assay procedure". All other kit reagents are ready-to-use. No reconstitution, mixing, or dilution is required.

**11. Assay procedure****Preparatory steps**

- (1) *Prepare an ethanol series (70%, 90%, and 100% ethanol solutions):* Dilute 100% ethanol with deionized or distilled water. These solutions can be stored in suitable containers and can be re-used.
- (2) *Heat Pretreatment Solution EDTA (PT2):* Heat to 98°C in a covered staining jar.
- (3) *Preparation of 1x Wash Buffer TBS:* Dilute 1 part of 20x Wash Buffer TBS (WB5) in 19 parts deionized or distilled water.  
*Diluted 1x Wash Buffer TBS is stable for one week when stored at 2-8°C.*
- (4) *1x Wash Buffer TBS:* For stringency wash, heat to 55°C in a covered staining jar.
- (5) *ZytoFast CISH Probe:* Before use, bring to hybridization temperature and mix thoroughly.
- (6) *Rabbit-Anti-DIG (AB11), Anti-Rabbit-AP-Polymer (AB12), Permanent Red Solution A (SB8a), Permanent Red Solution B (SB8b), Mayer's Hematoxylin Solution (CS1), Mounting Solution (alcoholic) (MT4):* Bring to RT (18°C-25°C) before use.

### Pretreatment (dewax/proteolysis)

- (1) Incubate slides for 10 min at 70°C (e.g., on a hot plate).
- (2) Incubate slides for 2x 5 min in xylene.
- (3) Incubate slides for 3x 3 min in 100% ethanol.
- (4) Wash slides 2x 1 min in deionized or distilled water at RT.
- (5) Incubate for 15 min in pre-warmed Heat Pretreatment Solution EDTA (PT2) at 98°C.

Use eight slides per staining jar (add dummy slides if needed).

- (6) Transfer slides immediately to deionized or distilled water and wash for 2x 2 min.
- (7) Apply (dropwise) Pepsin Solution (ES1) to the specimen and incubate for 5-15 min at 37°C in a humidity chamber.

*ES1 may form precipitates, which do not affect the quality.*

*As a general rule, we recommend to ascertain the optimum time for proteolysis in pre-tests.*

- (8) Immerse slides in deionized or distilled water at RT.
- (9) Dehydration in: 70%, 90%, and 100% ethanol, each for 1 min.
- (10) Air dry sections.

*Note: Make sure to completely dry sections prior to probe application.*

### Denaturation and hybridization

- (1) Pipette 10 µl of the ZytoFast CISH Probe onto each pretreated specimen.
- (2) Cover specimens with a 22 mm x 22 mm coverslip (avoid trapped bubbles) and seal the coverslip.

*We recommend using rubber cement (e.g., Fixogum) for sealing.*

- (3) Place slides on a hot plate or hybridizer and denature specimens for 5 min at 75°C.
- (4) Transfer slides to a humidity chamber and hybridize (e.g., in a hybridization oven) for 1 h at 37°C for DNA targeting\* probes or at 55°C for RNA targeting\* probes.

*\*Please refer to the package insert accompanying the probe. It is essential that specimens do not dry out during the hybridization step.*

### Post-hybridization and detection

- (1) Carefully remove the rubber cement or glue.
- (2) Remove the coverslip by submerging the slides in 1x Wash Buffer TBS at RT for 5 min.
- (3) Wash slides for 5 min in 1x Wash Buffer TBS at 55°C.

Use eight slides per staining jar (add dummy slides if needed).

- (4) Wash slides for 5 min in 1x Wash Buffer TBS at RT.
- (5) Apply Rabbit-Anti-DIG (AB11) (1-2 drops per slide) to the slides and incubate for 30 min at 37°C in a humidity chamber.
- (6) Wash slides 3x 1 min in 1x Wash Buffer TBS at RT.
- (7) Apply Anti-Rabbit-AP-Polymer (AB12) (1-2 drops per slide) to the slides and incubate for 30 min at 37°C in a humidity chamber.
- (8) Wash slides 3x 1 min in 1x Wash Buffer TBS at RT.
- (9) Prepare Permanent Red Solution (working solution): fill 1 ml Permanent Red Solution B (SB8b) in a graduated cup and add 16 µl Permanent Red Solution A (SB8a). Mix well.
- (10) Apply Permanent Red Solution (1-2 drops per slide) to the slides and incubate for 20 min at 37°C in a humidity chamber.
- (11) Wash slides 3x 1 min in deionized or distilled water at RT.
- (12) Counterstain specimens for 10-30 sec with Mayer's Hematoxylin Solution (CS1).
- (13) Transfer slides into a staining jar and wash 2 min under cold running tap water.
- (14) Dehydrate 3x 30 s in 100% ethanol (use very pure ethanol).
- (15) Incubate slides for 2x 30 s in xylene (use very pure xylene).
- (16) Air dry for approximately 2 min.
- (17) Avoiding trapped bubbles, cover the samples with a coverslip (22 mm x 22 mm; 24 mm x 32 mm) by using Mounting Solution (alcoholic) (MT4). Allow 20-30 min for the coverslip to become immobilized.

*Using a pipette tip which has been cut off to increase the size of the opening, can make the pipetting process easier.*

- (18) Evaluate stained specimens by light microscopy.

## 12. Interpretation of results

Using the ZytoFast PLUS CISH Implementation Kit AP-Permanent Red, hybridization of Digoxigenin-labeled oligonucleotides appear as bright red colored precipitates. Counterstaining the samples with Mayer's Hematoxylin Solution (CS1) will result in nuclei stained in blue.

Depending on the ZytoFast probe applied, a positive reactivity in the target cells, is found either within the cytoplasm or the nucleus, respectively. For a more detailed description of the expected signal pattern, please refer to the instructions of use accompanying the ZytoFast probe.

## 13. Recommended quality control procedures

Refer to the instructions for use of the respective ZytoVision probe.

## 14. Performance characteristics

Refer to the instructions for use of the respective ZytoVision probe.

## 15. Disposal

The disposal of reagents must be carried out in accordance with local regulations.

## 16. Troubleshooting

Any deviation from the operating instructions can lead to inferior staining results or to no staining at all.

### Weak signals or no signals at all

| Possible cause  | Action  |
|---|---|
| Cell or tissue sample has not been properly fixed   | Optimize fixing time and fixative   |
| Heat pretreatment, proteolysis, hybridization, denaturation, stringency wash or antibody-incubation temperature not correct | Check temperature of all technical devices used, using a calibrated thermometer. Use always the same number of slides in solutions with critical temperature  |
| Proteolytic pretreatment not carried out properly   | Depending on multiple factors, e.g., nature and duration of fixing, thickness of sections, and nature of tissue/cells, different incubation times may be required. Ascertain the optimum time for pepsin incubation in pre-tests  |
| Hybridization time too short  | Hybridize for at least 1 h; extend hybridization time if necessary  |
| Too low concentrated Wash Buffer  | Check concentration of Wash Buffer  |
| Old dehydration solutions   | Prepare fresh dehydration solutions   |
| Probe evaporation   | When using a hybridizer, the use of the wet stripes/water filled tanks is mandatory. When using a hybridization oven, the use of a humidity chamber is required. In addition, the coverslip should be sealed completely, e.g., with Fixogum, to prevent drying-out of the sample during hybridization |
| Insufficient preparation of chromogenic substrate   | check pipettes for proper performance; use calibrated pipettes  |
| Incubation temperature for color substrates not correct   | Check temperature of all technical devices used, using a calibrated thermometer   |
| Counterstaining time too long   | The counterstaining time depends on the nature of the specimen and should be optimized accordingly. Avoid dark counterstaining, because it may obscure positive staining signals  |

|  |   |
|--|---|
| No target sequences available                    | Use verified positive tissue to confirm test performance                                  |
| Bleuing of counterstain not carried out properly | Use cold running tap water for bleuing; do not use warm or hot water, or bleuing reagents |

**Signals too strong**

| Possible cause                                | Action   |
|---|--|
| Proteolytic pretreatment carried out too long | Depending on multiple factors, e.g., nature and duration of fixing, thickness of sections, and nature of tissue/cells, different incubation times may be required. Ascertain the optimum time for pepsin incubation in pre-tests |
| Substrate reaction is too intense             | Shorten substrate incubation time; do not heat substrate solution above temperatures given in the instructions for use   |

**Signals fade or merge**

| Possible cause                                | Action  |
|---|---|
| An unsuitable mounting solution has been used | Use only the mounting solution provided with the kit or xylene-based mounting solutions free of any impurities; do not use coverslip tape |

**Uneven or in some parts only very light staining**

| Possible cause   | Action  |
|--|---|
| Incomplete dewaxing  | Use fresh solutions; check duration of dewaxing times                   |
| Reagent volume too small                                   | Ensure that the reagent volume is large enough to cover the tissue area |
| Air bubbles caught before hybridization or during mounting | Avoid air bubbles   |

**Inconsistent results**

| Possible cause  | Action   |
|---|--|
| Insufficient drying before probe application  | Extend air-drying  |
| Too much water/wash buffer on tissue prior to application of pepsin, antibodies and/or color substrates | Ensure that excess liquid is removed from tissue section by blotting or shaking it off the slide. Small amounts of residual water/wash buffer do not interfere with the test |
| Variations in tissue fixation and embedding methods   | Optimize fixation and embedding methods  |
| Variations in tissue section thickness  | Optimize sectioning  |

**Tissue morphology degraded**

| Possible cause                                    | Action                            |
|---|-----------------------------------|
| Cell or tissue sample has not been properly fixed | Optimize fixing time and fixative |
| Proteolytic pretreatment not carried out properly | Optimize pepsin incubation time   |

**Noisy background**

| Possible cause                          | Action   |
|---|--|
| Stringency wash temperature not correct | Check temperature of the technical devices used, using a calibrated thermometer. Use always the same number of slides in the jar. We recommend not to use more than eight slides per jar for heat incubation steps |

|   |   |
|---|---|
| Slides not thoroughly rinsed                              | Use fresh and sufficient wash buffer and deionized or distilled water where indicated |
| Sections dried out any time during or after hybridization | Avoid sections being dried out; use humidity chamber; seal coverslip properly         |
| Prolonged substrate incubation time                       | Shorten substrate incubation time   |
| Incomplete dewaxing                                       | Use fresh solutions; check duration of dewaxing                                       |
| Proteolytic pretreatment too strong                       | Optimize pepsin incubation time   |
| Slides cooled to room temperature before hybridization    | Transfer the slides quickly to hybridization temperature                              |

**Overlapping signals**

| Possible cause                             | Action                                       |
|--|--|
| Inappropriate thickness of tissue sections | Prepare 3-5 $\mu\text{m}$ microtome sections |

**Specimen floats off the slide**

| Possible cause                      | Action                                      |
|-------------------------------------|---|
| Unsuitable slide coating            | Use appropriate (positively charged) slides |
| Proteolytic pretreatment too strong | Shorten pepsin incubation time              |

**17. Literature**

- Isola J, Tanner M (2004) *Methods Mol Med* **97**: 133-44.
- Speel EJ, et al. (1994) *J Histochem Cytochem* **42**: 1299-307.
- Tsukamoto T, et al. (1991) *Int J Dev Biol* **35**: 25-32.
- Wilkinson DG: *In Situ Hybridization, A Practical Approach*, Oxford University Press (1992), ISBN 0 19 963327 4.

Our experts are available to answer your questions. Please contact [helptech@zytovision.com](mailto:helptech@zytovision.com)



ZytoVision GmbH  
Fischkai 1  
27572 Bremerhaven/ Germany  
Phone: +49 471 4832-300  
Fax: +49 471 4832-509  
[www.zytovision.com](http://www.zytovision.com)  
Email: [info@zytovision.com](mailto:info@zytovision.com)

**Trademarks:**

ZytoVision® and ZytoFast® are trademarks of ZytoVision GmbH.