

VisionArray® Analysis Package SingleScan

For visualization and quick analysis of the VisionArray® Chips data

- Legend
- 1 Scanner 8100
 - 2 Hand Scanner
 - 3 Slide Holder
 - 4 Laptop
 - 5 USB-Hub
 - 6 External Hard Drive
 - 7 Computer Mouse



Prod. No. Product

E-4060-1

VisionArray Analysis Package SingleScan CE IVD

Incl. Scanner 8100; Slide Holder SingleScan; Hand Scanner; PC with preinstalled VisionArray Analyzer Software SingleScan; USB-Hub; External Hard Drive; Computer Mouse

VisionArray® Analyzer Software SingleScan

- Simple visualization and quick analysis of the VisionArray® Chip data
- Analysis of a Chip and the report of the results can be achieved in just a few minutes
- Program navigation is easy and intuitive for the user
- Scans are stored including all sample and Chip data in an integrated database on the enclosed external hard drive

VisionArray® Analyzer
19/09/2019 10:34:38

Start Archive x 0046_A001-SF01-001 x

Scan-ID: 46	Name:	Slide-ID: A001-SF01-001	Date of Scan: 19/09/2019 10:14
Case No.: 242	Date of Birth:	Chip-Lot: SF01	Chip File: VisionArray HPV Chip 1.0

Zoom: 12% | Overlay On/Off

Expert Mode:

✓ Hybridization/Grid
✓ Positive Control ✓ Duplicates

Target	Int. 1	Int. 2
HPV11	115.05	63.27%
HPV18	115.94	63.76%
HPV82 (IS39)	108.61	59.73%
Hyb./Grid	181.84	100.00%
Pos.	136.54	75.09%

Activity	Details	User	Time Sta
Scan created		Default	19/09/201
Automatic Analysis		Default	19/09/201

Next