

VisionArray® Analysis Package MultiScan

For visualization and quick analysis of up to 6 VisionArray® Chips simultaneously



Legend

- 1 Scanner V600 Photo
- 2 Laptop
- 3 USB-Hub
- 4 External Hard Drive
- 5 Computer Mouse

Prod. No.	Product
E-4070-1	VisionArray Analysis Package MultiScan CE IVD NEW Incl. Scanner V600 Photo; Slide Holder MultiScan; PC with preinstalled VisionArray Analyzer Software MultiScan; USB-Hub; External Hard Drive; Computer Mouse

VisionArray® Analyzer Software MultiScan

- Simple visualization and quick analysis of up to 6 VisionArray® Chips simultaneously
- All available VisionArray® Chips can be combined using the Scanner V600 Photo and are automatically detected by the software offering maximum flexibility
- Analysis of the Chips and the report of the results can be achieved in just a few minutes
- Scans are stored including all sample and Chip data in an integrated database on the enclosed external hard drive

CE IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.